

StayBright™ Plus

Professional Take-Home Whitening Enhancer



Why *StayBright™ Plus* Works

StayBright's formulation uses materials that have been in use in dentistry for over 100 years and are clinically proven to be safe and effective.

How *StayBright™ Plus* Works

Your patient simply brushes a thin layer of StayBright Take-Home Whitening Enhancer on each tooth to be whitened. StayBright incorporates the peroxide within a polymer base which leaves behind a film that temporarily adheres to the tooth surface. Proteins from the saliva on the tooth surface and inside the lips immediately act on the peroxide ions, causing an effervescent effect that releases the peroxide on the tooth surface. The peroxide then begins to whiten the teeth. The polymer film creates a barrier that prevents the saliva from breaking down the peroxide prematurely. The reaction lasts for approximately 20 minutes. The polymer compound can be easily brushed off, or it will wear off after an hour.

Safe, Quick and Affordable
The best choice for your patients' take-home whitening needs

2830 Via Orange Way, Ste B
Spring Valley, CA 91978
USA
619.660.5410 phone
619.660.5459 fax
800.400.2262 toll free
www.luma-lite.com
email: info@luma-lite.com

LUMALITE™

ISO 9001:2000
ISO 13485 Compliant
Certificate Number: A2386US

StayBright™ Plus

Professional Take-Home Whitening Enhancer
Item #3056

Easy To Use

- Your patient simply brushes the StayBright™ Plus Take-Home Whitening Enhancer onto their teeth, eliminating the need for mouth trays

No Trays, Less Hassle

- Typical tray-based systems require valuable chair and lab time. StayBright™ Plus provides a professional solution that eliminates the need for tray fabrication

No Tooth Sensitivity

- StayBright™ Plus has no reported sensitivity compared to other types of take-home whitening systems

Minimal Treatment Time (per application)

- Each application takes only seconds to apply, with 2-3 applications recommended per day

Sold in a case of 6
1/3 oz. bottles.

Includes a Point of Sale Display



Uses

- Enhance and maintain chair-side whitening
- Patients can whiten as many teeth as they like, when they like
- As easy to apply as nail polish
- Gives your patient the freedom to apply the whitening compound when convenient

StayBright™ Plus Take-Home Whitening Enhancer

Also provided in:

LumaWhite™ Plus and GentleBright™ Plus Single Patient Whitening Kits